

hurrimix



***AN EFFORTLESS
ACTIVITY***

Equipment
Automatic mixing unit

Zhermack[®]
CLINICAL



hurrimix

Zhermack, the market leader in solutions for the dental sector, is proud to present *hurrimix*, the **new, automatic mixer** for alginates and stones.

Hurrimix is one of the most modern automatic mixers currently available in the dental market, and produces **superb instant mixes**.

Compact and **silent** in operation, *hurrimix* guarantees **consistent mixing quality** and is the ideal choice for mixing *Zhermack* alginates and stones.



**Stylish
Design**

The presence of bubbles in mixes can cause distortion and other imperfections, compromising the quality of the finished work. *Hurrimix* ensures consistent homogenous mixes with minimal effort from the operator.

hurrimix

APPLICATIONS

Alginate mixing



Stone mixing*



*Guaranteed mixing quality for Zhermack's Elite® Master, Elite® Rock, Elite® Stone, Elite® Model, Elite® Model fast, Elite® Ortho, Elite® Arti stones. Maximum mixing quantity: 100g

CHARACTERISTICS



■ Satin steel body with coloured, soft touch plastic cover



■ Electronic microprocessor controller. Waterproof controls protected by a polycarbonate membrane



■ Safety mechanism ensures the lid has to be properly closed to prevent accidental opening and any loss of material during mixing



■ Max. container capacity: 80g for alginates, 100g for stones. The mixing containers have rounded internal surfaces to easily remove the mixed material and for cleaning

SOFTWARE



Extremely intuitive from the very first use. Large digital display for extremely simple, rapid operation.

20 mixing programs: 10 for alginates and 10 for stones, including 7 pre-set programs for type for Zhermack products. Users can modify any programs at any time.



Hurrimix boasts a compact, stylish design with harmonious lines that emphasise robustness and top quality materials.

hurrimix

CUSTOMIZATION

Customise your own *hurrimix* with 6 different coloured covers*.



* White (standard), orange, green, blue, navy, and grey (optional)



Hurrimix: a range of benefits in only a few seconds of effort

- Creamy mix
- Improved mechanical characteristics from alginates
- Quick, practical mixing
- Less material wasted
- Minimum user effort

SPOTLIGHT ON TECHNOLOGY





Technical specifications

Product code	Model	Dimensions	Weight	Mains voltage	Consumption
C305130	Hurrimix 230V, standard white cover	W24 / L30 / H31 cm	22 kg	230 Vac / 50 Hz	400 Watt
C305131	Hurrimix 115V, standard white cover	W24 / L30 / H31 cm	22 kg	115 Vac / 60 Hz	400 Watt

Standard equipment

Standard equipment		Accessories	
XR0250820	Mixing removal spatulas	XR0250860	Orange cover
XR0250830	Pre-mixing stick	XR0250865	Navy cover
XR0250800	Alginate mixing containers x2	XR0250870	Green cover
XR0250810	Zhermack stone mixing containers x2	XR0250855	Blue cover
XR0250840	Water bottle with graduated measuring cup	XR0250850	White cover
		XR0250875	Grey cover



For further information on this product, visit the Zhermack website www.zhermack.com

BRANCHES

GERMANY
Zhermack GmbH
Deutschland
Öhlmühle 10 - D-49448 Marl
Tel. +49 - 0 54 43 / 20 33 - 0
Fax +49 - 0 54 43 / 20 33 - 11
info@zhermack.de
www.zhermack.com

POLAND
Zhermapol Sp. z o.o.
ul. Augustówka 14
02 - 981 Warszawa - Polska
Tel. +48 - 22 858 82 72
+48 - 22 858 73 41
Fax +48 - 22 642 07 14
biuro@zhermapol.pl
www.zhermapol.pl

USA
Zhermack Inc.
P.O. Box 4195, River Edge
New Jersey 07661 - 4195
Tel. +1 - (877) 819 6206
+1 - (732) 389 8540
Fax +1 - (732) 389 8543
info@zhermackusa.com
www.zhermackusa.com

COUNTRY CONTACTS

UK AND IRELAND
Tel. +44 - (0)7870 690 811
uk@zhermack.com

SPAIN
Tel. +34 - 900 99 39 52
info.spain@zhermack.com

FRANCE
Tel. +33 - 0800 915 083
info.france@zhermack.com

RUSSIA
Tel. +7 - 916 930 1191
info.russia@zhermack.com

MEXICO AND CENTRAL AMERICA
Tel. +52 - 1 55 919 804 97
info.mexico@zhermack.com

SOUTH AMERICA
Tel. +57 - 310 2601 777
info.colombia@zhermack.com

Zhermack SpA

Via Bovazecchino, 100 - 45021 Badia Polesine (RO) ITALY Tel. +39 - 0425 597611 Fax +39 - 0425 53596 info@zhermack.com www.zhermack.com