

For further information on these products please visit the website en.zhermack.com



elite[®] LC tray **ROUND**

DESIGNED FOR EFFICIENCY



Light Curing Tray Material

Fruit of Zhermack's decade of know-how in making light curing plates, here for you is the "ROUND" version of the **elite[®] LC tray** line.

elite[®] LC tray ROUND plates have an **innovative ergonomic shape**, designed to reduce the amount of material left after working to the minimum. The prosthodontist no longer has to worry about leftover material, saving time and avoiding possible waste.

elite[®] LC tray ROUND is therefore the ideal product for laboratories wanting medium-sized ergonomically shaped plates, **optimised to improve efficiency**.



en.zhermack.com
info@zhermack.com

Zhermack[®]
BEYOND INNOVATION

Benefits

- Ergonomic shape to reduce leftover material
- Imperceptible contraction after curing
- Easy to work, does not stick to the hands
- Hardness and rigidity designed to facilitate finishing
- Ideal thickness for multiple applications
- Shelf-life 3 years

USES



Individual and functional impression trays



Bite plates



Bases for mounting teeth and/or fixing registration devices



Blocking telescopic or secondary crowns

Technical data

Curing times	3' - 5'
Thickness	2,4 mm*
Hardness at 1' (Shore D)	80 Shore D
Curing	UV (350-400 nm) or halogen (420-480 nm) light

* tolerance +/- 0,1 mm

Packs

D500050	Elite® LC Tray ROUND White (50 pcs)
D500052	Elite® LC Tray ROUND Pink (50 pcs)
D500054	Elite® LC Tray ROUND Blue (50 pcs)

BRANCHES

GERMANY
Zhermack GmbH
Deutschland
Öhlmühle 10 - D-49448 Marl
Tel. +49 - 0 54 43 / 20 33 - 0
Fax +49 - 0 54 43 / 20 33 - 11
info@zhermack.de
www.zhermack.com

POLAND
Zhermapol Sp. z o.o.
ul. Augustówka 14
02 - 981 Warszawa - Polska
Tel. +48 - 22 858 82 72
+48 - 22 858 73 41
Fax +48 - 22 642 07 14
biuro@zhermapol.pl
www.zhermapol.pl

USA
Zhermack Inc.
P.O. Box 4195, River Edge
New Jersey 07661 - 4195
Tel. +1 - (877) 819 6206
+1 - (732) 389 8540
Fax +1 - (732) 389 8543
info@zhermackusa.com
www.zhermackusa.com

COUNTRY CONTACTS

UK AND IRELAND
Tel. +44 - (0)7860 346 129
uk@zhermack.com

SPAIN
Tel. +34 - 900 99 39 52
info.spain@zhermack.com

TURKEY
Tel. +90 533 76 44 955
info.turkey@zhermack.com

FRANCE
Tel. +33 - 0800 915 083
info.france@zhermack.com

RUSSIA
Tel. +7 - 916 930 1191
info.russia@zhermack.com

MEXICO AND CENTRAL AMERICA
Tel. +52 - 1 55 919 804 97
info.mexico@zhermack.com

SOUTH AMERICA
Tel. +57 - 310 2601 777
info.colombia@zhermack.com

Zhermack SpA

Via Bovazecchino, 100 - 45021 Badia Polesine (RO) ITALY
Tel. +39 - 0425 597611 Fax +39 - 0425 53596 info@zhermack.com www.zhermack.com