

Elite Double. Professional solutions for model duplications

Elite Double, the newly improved range of Zhermack duplication silicones, satisfies all dental laboratory needs.

Thanks to re-styling of the entire range, Elite Double is proud to offer easy handling and the ability to supply continuously uniform results in all situations. Its 1:1 ratio and snap set effect permit dental technicians to create precise, fast duplications with the added advantage of maintaining adequate working times.

This new range is composed of four specific hardness and three distinct setting times:

Elite Double 8 NORMAL – for core duplications in whole ceramic coatings

Elite Double 16 FAST / EXTRA FAST – for core duplications in whole ceramic coatings, for plaster and/or resin model duplications and for model duplications in controlled expansion coatings

Elite Double 22 NORMAL / FAST / EXTRA FAST – for plaster and/or resin model duplications and for model duplications in controlled expansion coatings

Elite Double 32 NORMAL / FAST / EXTRA FAST – for plaster and/or resin model duplications and for model duplications in free expansion coatings

In addition to previous ranges, Zhermack has introduced the new Elite Double 16 in FAST and EXTRA FAST versions, Elite Double 22 and 32 in EXTRA FAST versions and new 10kg packs, where available.

For the first time on the market, Zhermack presents the **EXTRA FAST** version, formulated especially for the **Doublemix** automatic mixer, an upcoming innovative mixing system for duplication and model removal in just 6 minutes.

Zhermack SpA

Via Bovazecchino 100
45021 Badia Polesine RO - Italy
Tel. 0425 597611
Fax 0425 53596
e-mail: info@zhermack.com
www.zhermack.com